




LOW SPEED CENTRIFUGE

LC500



TECHNICAL DATA



Program memory



Temperature protection & overshock protection



Multi rotors optional



Double electronic lock



Model	LC500-8	LC500-12	LC500-6	LC500-18	LC500-24
Rotor Capacity	15ml*8 Centrifuge tube 2/5/10ml*8 Blood tube	15ml*12 Centrifuge tube 2/5/10ml*12 Blood tube	15/50ml*6 Centrifuge tube 2/5/10ml*6 Blood tube	2/5/10ml*18 Blood tube	2/5/10ml*24 Blood tube
Speed	500 - 4000rpm				
Max RCF	1967xg	2075xg	2146xg	2182xg	2164xg
Semidiameter Of Rotor	11cm	11.6cm	12cm	12.2cm	12.1cm
Timing Range	30s - 59min50s				
Protection Function	Lid opening protection, cavity overtemperature protection, vibration projection				
Power	60W	60W	80W	80W	80W
Power Supply	AC 110V - 220V, 50/60Hz				



**Standard
MEDICAL**

Av Anatole France Imm Intilaka Centre A1-1, Sfax, Tunisia.
Tél.: 58 520 400, 58 520 402, 58 520 404.
E-mail : standardmedical7@gmail.com.

