

High Precision Constant Temperature Incubator

© High Precision bacterial microbial incubator.



DH45L, DH65L, DH125L, DH210L

Product Description

- 1.** The double door structure adopts high quality glass, easy to observe samples, the outdoor with magnetic strip is easy to be opened and closed, good sealing, open door power-off function.
- 2.** Standard lights, UV lamps, with breeze circulation fan.
- 3.** High precision digital independent temperature limiter.
- 4.** With function of cutting off UV lamp by open door.



DH45L, DH65L, DH125L, DH210L

Model		DH45L	DH65L	DH125L	DH210L
Cycle Mode		Breeze circulation			
Function	Tem. Range	AT+5-70°C			
	Term. Resolution Rato	0.1°C			
	Tom. Motion	±0.5°C			
	Tem. Uniformity	+08 °C			
Structure	Inner Chamber	Mirror stainless steel			
	Outer Shell	Cold rolling steel electrostatic spraying exterior			
	Insulation Layer	High quality foam board			
	Heater	Mica electrothermal fim			
	Power Rating	0.35KW	0.45KW	0.6kW	0.7kW
	Exhaust Hole	Ø28mm top with function of test hole)			
Controller	Tem. Control Mode	PID Intelligent			
	Tem. Setting Mode	Touch button soting			
	Tem. Display Mode	Digital temperature: LCD upper row, Setting temperature: the lower row			
	Timer	0-9999 min (with timing wait function)			
	Operation Function	Fixed temperature operation, timing function, auto stop.			
	Additional Function	Breeze Circulation fan, UV lamp,LED Floodlight, Sensor deviation correction, temperature overshoot self-tuning power-off parameter memory			
	Sensor	PT100			
Safety Device		Digital, over temperature sound-light alarm			
Specification	Inner Chamber Size(W"L"HXmm)	350 350 350	400*350*450	500*450*550	600'580'600
	Exterior Size (W"L"H)(mm)	525'480'620	575°480*720	675°580*820	775'710'870
	Packing Size (W*L*H)(mm)	605'572'775	656'572'876	775'710'870	855 802*1025
	Volume	45L	65L	125L	210L
	Shelf Number	7	9	13	14
	Load Per Rack	15kg			
	Shelf Space	35mm			
	Supply(50/60HZ)Current Rating	AC220V/1.6A	AC220V/2.0A	AC220V/2.7A	AC220V/3.2A
	NW/GW (kg)	27/30	32/35	45/49	58/63
	Shelf	2			
	Shelf Frame	4			
Optional Accessories		Shelf, U-Disk storage, Printer, Remote Control, Data Recorder			



**Standard
MEDICAL**

Av Anatole France Imm Intilaka Centre A1-1, Sfax, Tunisia.
Tél.: 58 520 400, 58 520 402, 58 520 404.
E-mail : standardmedical7@gmail.com.

